



SEIDENADER
KÖRBER SOLUTIONS

THE WHOLE FORCE IN SMALL SIZE.



HSA O₂


Leak testing for lyo products on smallest footprint. The battle against invisible leakages continues.

HSA O₂ – YOUR KEY BENEFITS

- ▶ Easily integratable into standard inspection machines
 - no extra space required for a separate machine
 - lower investment costs
- ▶ Only one sensor needed. This guarantees lower validation effort and reproducibility
- ▶ High speed analysis up to 24,000 containers per hour

RELIABLE INSPECTION RESULTS FOR

- ▶ Vials with ≥ 16 mm inside diameter (2 R/10 R/15 R)
- ▶ Filled and sealed under protective atmosphere



Leak diameter
in μm

Waiting time in h until balance
through diffusion (21% O₂)

5

300

10

90

20

25